


**Experiment 25: Determination of  $K_c$  for a Chemical Reaction**

*This experiment was adapted from Ross S. Nord, Chemistry Department, Eastern Michigan University.*

**Purpose:** To determine the equilibrium constant for the reaction between iron(III) and thiocyanate, and to study the equilibrium system of this chemical reaction.

**Background:** Iron(III) and thiocyanate react to form an equilibrium mixture with the product formed; the iron-thiocyanate complex ion. The net-ionic equation is:



The iron(III) solution is slightly pale-yellow in color, the thiocyanate solution is colorless, and the complex ion product is an orange-red color. The colored product will allow for the use of a spectrophotometer to quantitate the amount of product made. Using the data from the spectrophotometer, along with concentrations of iron(III) and thiocyanate, we will study this equilibrium system.

When a reaction reaches equilibrium, not all of the reactants became product. A mixture of reactants and products exist. The equilibrium constant is calculated as the ratio of the concentration of product(s) divided by the concentration of reactant(s). This is the simplified version of the math. The more complicated version involves calculations with activity coefficients and ionic strength due to the fact that this equilibrium system involves ionic species. Overall, we can say that it is typical for the math to be simplified for the calculation of the equilibrium constant.

$$K_c = \frac{[\text{FeNCS}^{2+}]}{([\text{Fe}^{3+}] \times [\text{SCN}^{1-}])}$$

This experiment consists of three parts:

**Part 1:** Determine the  $\lambda_{\text{max}}$  for the  $\text{FeNCS}^{2+}$  ion. Since the solution is an orange-reddish color, wavelengths in the 400 nm to 500 nm range are being absorbed.

**Part 2:** The reaction between the iron(III) and thiocyanate will be studied as an equilibrium, with the ion concentration as low as possible. Reasonable amounts of each reactant will be combined to allow the equilibrium mixture to form.

**Part 3:** The reaction between the iron(III) and thiocyanate will be studied as an equilibrium, with the ion concentration increased by the addition of 0.10 M  $\text{KNO}_3$ . The amount of each reactant will be the same as in Part 2.

**Chemicals** *(Use caution with each of these solutions.)*

0.0020 M Fe(NO<sub>3</sub>)<sub>3</sub> in 0.0050 M HNO<sub>3</sub>

0.0020 M NaSCN in DI water

0.10 M KNO<sub>3</sub>

Deionized water

**Equipment**

1.0 mL pipet (1), 10.0 mL pipet (2), pipettor (2), beakers (3)

Small test tubes

Parafilm

Spectrophotometer

**Procedure**

**Part 1: Determining  $\lambda_{\max}$  for FeNCS<sup>2+</sup>**

- a) Combine 1.50 mL of Fe<sup>3+</sup> with 2.00 mL of SCN<sup>1-</sup> in a small test tube.
- b) Cover the test tube with Parafilm. Secure your thumb over the Parafilm and invert the test tube to mix three times. This is the equilibrium mixture.
- c) Set the software on the spectrophotometer to perform a wavelength scan.
- d) Place a test tube of Fe<sup>3+</sup> solution in the sample compartment of the spectrophotometer as the blank and zero the.
- e) Place the test tube with the equilibrium mixture in the sample compartment and run the wavelength scan.
- f) The spectrophotometer will display a graph of wavelength vs. absorbance. Record the wavelength of maximum absorbance.

**Part 2: Observe the Equilibrium Mixture with Low Ionic Strength**

- a) Prepare a series of equilibrium mixtures, each in its own small test tube. Use the volumes of each reactant listed in Table 1.
- b) Cover each test tube with Parafilm. Secure the film with your thumb and invert to mix three times.
- c) Set the spectrophotometer wavelength to  $\lambda_{\max}$ . Use the Fe<sup>3+</sup> solution as your blank and zero the spectrophotometer.
- d) Measure the absorbance of each equilibrium mixture.
- e) Add 1.5 mL of DI water to each test tube and blank, and invert to mix three times.
- f) Rezero with the blank and measure the absorbance of each equilibrium mixture.
- g) Add 1.0 mL of DI water to each test tube and blank, and invert to mix three times.
- h) Rezero with the blank and measure the absorbance of each equilibrium mixture.

**Table 1 Solution Preparation for Part 2**

Solution Number	Volume of Fe <sup>3+</sup> , mL	Volume of SCN <sup>1-</sup> , mL
Blank	4	0
1	0.50	3.00
2	1.00	2.50
3	1.50	2.00
4	2.00	1.50
5	2.50	1.00
6	3.00	0.50

**Part 3: Observe the Equilibrium Mixture with High Ionic Strength**

- Prepare a series of equilibrium mixtures, each in its own small test tube. Use the volumes of each reactant listed in Table 1 (same as in Part 1).
- Cover each test tube with Parafilm. Secure the film with your thumb and invert to mix three times.
- Set the spectrophotometer wavelength to  $\lambda_{\max}$ . Use the Fe<sup>3+</sup> solution as your blank and zero the spectrophotometer.
- Measure the absorbance of each equilibrium mixture.
- Add 1.5 mL of 0.10 M KNO<sub>3</sub> to each test tube and blank, and invert to mix three times.
- Rezero with the blank and measure the absorbance of each equilibrium mixture.
- Add 1.0 mL of 0.10 M KNO<sub>3</sub> to each test tube and blank, and invert to mix three times.
- Rezero with the blank and measure the absorbance of each equilibrium mixture.

**Calculations**

Use Excel or Google Sheets for the calculations. This will save a lot of time and provide experience with programming spreadsheets. There is a video in your Blackboard shell that shows you how to make the Excel / Google Sheets file for this experiment.

**You will have to calculate the following:**

- The initial concentration of iron(III).
- The initial concentration of thiocyanate.
- The concentration of FeNCS<sup>2+</sup> made (this is the equilibrium concentration). Remember Beer's Law,  $A=abc$ . Make the spreadsheet solve for c, the concentration. You know the absorbance, A. Use 5800 M<sup>-1</sup>·cm<sup>-1</sup> for the extinction coefficient a, and use 1 cm for the path length b.
- The equilibrium concentration of iron(III).
- The equilibrium concentration of thiocyanate.
- The value of the equilibrium constant, using equilibrium concentrations.

*All of these calculations are shown in the prelab PowerPoint file and in the video.*

### **Conclusion**

When you write your conclusion, consider these topics:

- a) Was the equilibrium constant a constant within each table, considering experimental error?
- b) Was the equilibrium constant a constant after the dilution of water, considering experimental error?
- c) Was the equilibrium constant affected by the addition of  $\text{KNO}_3$ , more so than the addition of DI water?
- d) The preparation of the original solutions was the same in Part 2 and Part 3. How do those equilibrium constant values compare?
- e) Did you notice any patterns with the data?