


**Experiment 26: Titration of a Commercial Antacid***(This experiment was adapted from CUNY.)***Purpose**To determine the amount of  $\text{CaCO}_3$  present in a commercial antacid tablet.**Background**

The parietal cells in the stomach secrete hydrochloric acid (HCl) at a concentration of roughly 0.16 M. The flow of HCl increases when food enters the stomach. If you eat or drink too much, you may develop heartburn or indigestion. Antacids are used to neutralize this excess acid.

The active ingredient in antacids is calcium carbonate,  $\text{CaCO}_3$ , a base. There are also other ingredients, such as binders and flavors present in each tablet.

HCl is neutralized by calcium carbonate as illustrated below:



To determine the quantity of  $\text{CaCO}_3$  in the antacid tablet, we are first going to dissolve the tablet in an excess amount of acid, HCl, of known concentration. Some of the HCl will be neutralized by the carbonate, but there will be some remaining. We will then perform a titration with NaOH to figure out the amount of excess HCl. Then, from this, we can calculate how much acid reacted with the antacid, followed by the amount of antacid present in the mass of tablet used. This method of analysis is called a back-titration.

The reactions above are reversible, which means that the  $\text{CO}_2$  dissolved in water would produce some carbonic acid. This acid would react with the NaOH we are titrating and give us inaccurate results. Therefore, it is important to boil the solution when the calcium carbonate reacts with acid, to remove  $\text{CO}_2$  as a gas. The net result is below:

**Chemicals**

0.10 M hydrochloric acid, HCl

0.10 M sodium hydroxide, NaOH

Antacid tablets, 500 mg  $\text{CaCO}_3$ 

Phenolphthalein, 1%

Deionized water

**Equipment**

Buret, buret stand and clamp

Erlenmeyer flask, 125 mL

Graduated cylinder, 50 mL or 100 mL

Hot plate (use the one in the hood)

Mortar &amp; Pestle

Beakers, as needed

**Procedure** (*Do one trial at a time, start to finish.*)

1. Rinse your buret with DI H<sub>2</sub>O, and then rinse it three times with 2 mL-portions of the NaOH solution. Fill the buret with the 0.10 M NaOH solution. Allow enough of the NaOH solution to drain from the buret to remove the air bubble in the tip of the buret. Refill the buret as needed.
2. There will be hot-water baths set up in the fume hoods. Let the bath water get hot enough to boil gently. You will be using this to heat your acid and antacid mixture.
3. Obtain two antacid tablets. Find their exact, combined mass with the lab balance.
4. Crush both tablets with a mortar and pestle, until they are a fine powder.
5. Use the graduated cylinder to add 50.0 mL of the 0.10 M HCl to the 125 mL Erlenmeyer flask. Rinse the graduated cylinder with 10 mL of DI water and add this rinse to the Erlenmeyer flask.
6. Determine the exact mass of approximately 0.450 g of this powder. Your sample mass can be slightly less than 0.450 g, but do not use more than 0.450 g. Transfer this sample into the Erlenmeyer flask.
7. Swirl the flask to dissolve the antacid powder in the acid and place this flask in the hot-water bath for 10 minutes. After the 10 minutes, remove the Erlenmeyer flask from the hot-water bath and allow it to cool on your lab bench. Placing the Erlenmeyer in a cool-water bath will help cool it faster.
8. Add three drops of phenolphthalein indicator to the Erlenmeyer flask. Mix well.
9. Record the initial buret volume. Titrate the solution in the flask until the endpoint is reached. The endpoint will be visible when the indicator turns very faint pink. Record the final buret volume.
10. Repeat steps 5 – 9 for two more trials. You should have a total of 3 useable trials. Do one trial at a time, from start to finish.

**Data Needed**

Exact mass of two antacid tablets

For each trial, make sure you record the following information:

- a) Volume of HCl added into the Erlenmeyer flask
- b) Exact mass of antacid powder added to the Erlenmeyer flask
- c) The initial volume of NaOH (before the titration is started)
- d) The final volume of NaOH (after the titration is stopped)

**Calculations**

For each trial, calculate the following: (remember,  $M = \text{moles/L}$  so  $M \times L = \text{moles}$ )

- Calculate the number of moles of HCl added into the Erlenmeyer flask
- Calculate the number of moles of NaOH delivered from the buret
- Subtract the moles NaOH from the moles HCl; this difference is the HCl used to neutralize the  $\text{CaCO}_3$ ; remember ( $2 \text{ H}^+ : 1 \text{ CO}_3^{2-}$ )
- Convert the moles of  $\text{CaCO}_3$  in the flask to grams of  $\text{CaCO}_3$
- Divide the grams of  $\text{CaCO}_3$  by the mass of the antacid powder put into the flask;  
( $\text{g CaCO}_3 / \text{g antacid powder in flask}$ )
- Divide the mass of the two antacid tablets by two; this gives the g antacid powder per one tablet
- Calculate the g  $\text{CaCO}_3 / \text{tablet}$  with the following calculation:

$$\frac{\text{g CaCO}_3}{\text{g antacid powder in flask}} \times \frac{\text{g antacid powder}}{\text{one tablet}} = \frac{\text{g CaCO}_3}{\text{tablet}}$$

- Average the three values of the g  $\text{CaCO}_3 / \text{tablet}$
- How does this average g/tablet compare to the manufacture's claim? Calculate the percent difference between your experimental average result and the value on the bottle of antacid. Use the average experimental result for this calculation.

$$\% \text{ difference} = \frac{/ \text{ true g/tablet} - \text{exp g/tablet} /}{\text{true g/tablet}} \times 100$$

\*the true g/tablet is the value from the bottle